signum [®] composite / composite flow / matrix [®]	
	
DE Gebrauchsanweisung	SE Bruksanvisning
GB Instructions for use	NO Bruksanvisning
FR Mode d’emploi	GR Οδηγίες χρήσης
ES Instrucciones de uso	CZ Návod k použití
IT Istruzioni per l’uso	HU Használati utasítás
PT Instruções de uso	LV Lietošanas instrukcija
NL Gebruiksaanwijzing	LT Naudojimo instrukcija
PL Instrukcja obsługi	
HR Upute za uporabu	
RO Instrucțiuni de utilizare	
KR 사용 설명서	
BG Инструкции за употреба	

Signum® composite / composite flow / matrix®

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Signum® composite / composite flow / matrix®

The products of the brand Signum are light curing veneering composites for crowns and bridges.

Fields of Application

- Full and partial veneering for fixed metal-supported crowns and bridges.
- Veneering of resin based frames (PMMA).
- Veneering of removable prostheses (telescope/coronic crowns and attachments).
- For the veneering of Implant supra-/tertiary structures.
- Zirconia-supported longterm provisionals.
- Temporary anterior, posterior crowns (w/o frame).
- Characterisation of artificial teeth.

Note

The metal supported veneering is permitted up to the second premolar teeth without limitation, unless the molars are expected to be subject to undue stress. Signum composite, Signum composite flow and Signum matrix may also be applied for the occlusal veneering of molars to obtain aesthetically pleasing results. In case of any foreseeable physiological strain or forces acting on the masticatory surfaces of the molars (e.g. bruxism or parafunctions), veneer contact must be ruled out or metal stops applied. Fully veneered pontics should be separated interdentially prior or after the polymerisation in order to minimise the risk of stress cracks. Especially for the press to – or injection technique a sufficient interdental separation for long-span bridges must exist!

Composition

Signum composite
Signum dentine, Signum margin, Signum enamel, Signum value, Signum effect: Multifunctional methacrylic acid esters, photo initiators, stabilising agents, inorganic pigments (cadmium free).

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