# signum<sup>®</sup>

## insulating gel

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GB GB Instructions for use

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## GB Instructions for use

Gel for insulation of composites and avoidance of an inhibition layer.

### Composition

Contains glycerine and silicic acid

### Area of application

Signum insulating gel is for insulation of plaster against composite materials and for avoidance of an inhibition layer during polymerisation in a laboratory light unit.

### Side effects

None known

### Reciprocal effects with other materials

None known

### Dosing and application

#### 1. Insulation

Two layers of Signum insulating gel should be brushed carefully onto the plaster to be insulated. Thereafter place the composite onto the insulated surfaces and cure with a light unit. To guarantee risk-free removal of the composite workpiece from the plaster model, undercuts within the preparation must be blocked out. Plaster models blocked out with resin require a special insulation. Signum insulating gel prevents adhesion. It is therefore imperative to remove all Signum insulating gel residues carefully with water prior to cementation of the workpiece.

#### 2. Avoidance of an inhibition layer

Signum insulating gel prevents the formation of a dispersion layer. A thin, incomplete polymerised layer, the dispersion layer, is formed on the surface of all composites under the action of atmospheric oxygen. This dispersion layer is a pre-requisite for chemical bonding between the individual components while layers are being built up.

After working and polishing, the surface of the veneer must be free of a dispersion layer. When Signum insulating gel is used, the surface also cures perfectly. No dispersion layer is formed. This considerably facilitates the finishing, which presents problems particularly in regions that are difficult to ac-cess. The completely polymerised surface affords additional protection as far as the prevention of plaque deposits and colour changes are concerned.

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### Procedure

– Layer and polymerise the composite in the usual manner, being careful to adhere to the stipulated polymerisation times.

– Apply a layer of gel about 0,5 mm thick over the entire veneer before endpolymerisation.

– The time for the final polymerisation is 180 sec.

– The gel does not cure during the polymerisation.

– Rinse the gel off with water before finishing.

### Note

To guarantee perfect polymerisation of the paste beneath the gel, the gel must not be contaminated. The container should therefore be closed tightly immediately after use.

### Storage

Do not store above 25°C (77°F). Do not use after expiry date. Avoid direct exposure to sun. Always keep container closed.

### Disposal instructions

Observe the safety data sheet or national regulations for disposal.

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® = registered trademark of Kulzer GmbH
® = marque déposée de Kulzer GmbH